

| Company product code | Wound Care Product Name | Qty/Box |
|--|--|---------|
| Exudate Management: Foam, Silicone: Absorption of mod-Irg amt of exudate & skin tears (silicone foam). Do not use with hydrogels. | | |
| Note: Foams are not covered for use for protection unless exception is provided; do not use Foam on feet. | | |
| Molylycke | | |
| 581011 | Mepilex Border Lite 4cm x 5cm | 10/box |
| 581200 | Mepilex Border Lite 7.5cm x 7.5cm | 5/box |
| 581300 | Mepilex Border Lite 10cm x 10cm | 5/box |
| 595011 | Mepilex Border Flex 12.5cm x 12.5cm | 10/box |
| 595311 | Mepilex Border Flex 10cm x 10cm | 10/box |
| 595411 | Mepilex Border Flex 15cm x 15cm | 10/box |
| 595611 | Mepilex Border Flex 15cm x 20cm | 10/box |
| 284100 | Mepilex Lite Non-Border 10cm x 10cm | 5/box |
| 294100 | Mepilex Non-Border 10cm x 10cm | 5/box |
| 294300 | Mepilex Lite Non-Border 15cm x 15cm | 5/box |
| Foam Post-Op Silicone dressing Use: Post-op wound | | |
| 496400 | Mepilex post op Border 10cm x 20cm | 10/box |
| 496600 | Mepilex post op Border 10cm x 30cm | 10/box |
| Gelling Fiber Use: Packing of wounds and to provide additional absorption under a foam dressing in mod-heavily exudating wounds. | | |
| Convatec | | |
| 420127 | Aquacel Hydrofiber Ribbon 1cm x 45cm | 5/box |
| 420672 | Aquacel Hydrofiber Extra 10cm x 10 cm | 10/box |
| Alginate Use: Packing of wound and to provide additional absorption under a foam dressing in mod exudating wounds | | |
| Smith & Nephew | | |
| SN66000519 | Algisite M Calcium Alginate Dressing 5cm x 5cm | 10/box |
| SN66000520 | Algisite M Calcium Alginate Dressing 10cm x 10cm | 10/box |
| SN66000521 | Algisite M Calcium Alginate Dressing 15cm x 20cm | 10/box |
| SN66000522 | Algisite M Calcium Alginate Dressing 2cm x 30cm | 5/box |
| Antimicrobial (10% Povidone Iodine non-adherent sheet) Use: For wounds with minimal exudate or to assist in the drying of wounds | | |
| 3M | | |
| P01481 | Inadine Dressing 5cm x 5cm | 25/box |
| P01512 | Inadine Dressing 9.5cm x 9.5cm | 25/box |
| Cadexomer Iodine Use: For sloughy moist wounds, effective up to 3 days | | |
| Smith & Nephew | | |
| SN66060630 | Iodosorb Ointment 10gm tube | 4/box |
| SN66001286 | Iodosorb Powder 3gm packet | 7/box |

| | | |
|--|--|----------|
| Gelling Hydrofiber Use: For mod-lrg amt of exudate. Caution: Do not use on dry wounds or wounds with minimal wound exudate. | | |
| Convatec | | |
| 413571 | Aquacel Ag+ Ribbon 2cm x 45cm | 5/box |
| 413567 | Aquacel Ag+ Extra 10cm x 10cm | 10/box |
| 0.2% Polyhexamethylene Biquanide (PHMB)) Use: Gauze dressing for wound packing including tunnels; open surgical wounds | | |
| Cardinal | | |
| Z7831AMD | AMD PHMB Packing Strips 0.63cm x 0.91.4m | 10/box |
| Z7832AMD | AMD PHMB Packing Strip 1.27cm x 0.91.4m | 10/box |
| Z7833AMD | AMD PHMB Packing Strip 2.51cm x 0.91.4m | 10/box |
| Z3332 | AMD Kerlix Rolls 11.4cm x 3.7m | 100/box |
| Charcoal and Silver Use: For malodourous wounds contains Silver. NOTE: do not cut | | |
| 3M | | |
| MAS105 | Actisorb Ag 220 Charcoal 10.5cm x 10.5cm | 10/box |
| Protective Cover Dressings: Use: Cover and protect areas. Note: Does not promote moist wound healing | | |
| Smith & Nephew | | |
| 66000317 | Primapore 8cm x 10cm | 20/box |
| 66000318 | Primapore 8cm x 15cm | 20/box |
| 66000319 | Primapore 10cm x 20cm | 20/box |
| 66000320 | Primapore 10cm x 25cm | 20/box |
| 660007135 | Primapore 6cm x 8cm | 50/box |
| 660007140 | Primapore 10cm x 35cm | 20/box |
| Superabsorbent (cover dressing) Use: For highly exuding wounds, has extra fluid handling absorbency capacity | | |
| 3M | | |
| PRD-500-050 | Kerramax 10cm x 10cm | 10/box |
| PRD-500-100 | Kerramax 13.5cm x 15.5cm | 10/box |
| PRD-500-120 | Kerramax 10cm x 22cm | 10/box |
| PRD-500-240 | Kerramax 20cm x 22cm | 10/box |
| PRD-500-380-B10 | Kerramax 20cm x 30cm | 10/box |
| Skin Protection: Use: Skin Barrier Wipes, Adhesive Remover, Wound filler Zinc Paste, Silver Transfer Sheet | | |
| Smith and Nephew | | |
| 66800709 | Skin Barrier Spray - no sting Skin Prep (28 ml bottle) | 1 bottle |
| 59420600 | Skin Barrier Wipes - no sting Skin Prep (1ml) | 50/box |
| 3M | | |
| 3M5050 | Cavilon Advanced (2.7ml applicator) | 20/case |
| Coloplast | | |
| COL1964 | Triad Hydrophillic Wound Paste, 71G | 12/box |
| COL7912 | Interdry Silver Roll 10 inches x 36 inches | 10/box |

| | | |
|---|---|--------|
| COL7910 | Interdry Silver Roll 10 inches x 144 inches | 1 Roll |
| Compression Use: For venous or mixed disease leg ulcers - suggest Drs orders and lower limb assessment before application. Once healed, compression stockings are suggested for long term maintenance. | | |
| Medical Mart | | |
| CA600001 | Edema Wear Small up to 45cm/ Navy (foot to knee) | 2/pk |
| CA960001 | Edema Wear Medium 45- 75cm/Yellow (foot to knee/groin) | 1/pk |
| Non-Adherent Contact Layer Petrolatum & Silicone Use: To prevent wound adherence | | |
| 3M | | |
| 2012 | Adaptic Petrolatum Impregnated Dressing 7.5cm x 7.5cm | 24/box |
| 2013 | Adaptic Petrolatum Impregnated Dressing 20cm x 7.5cm | 24/box |
| TCH502 | Adaptic Petrolatum Impregnated Dressing 11cm x 7.6cm | 10/box |
| TCH503 | Adaptic Petrolatum Impregnated Dressing 12.7cm x 15cm | 10/box |
| Non-Adherent (Digit-silicone) Use: Dressing designed for toes and digits | | |
| 3M | | |
| MAD042 | Adaptic Digit Non-Adherent Dressing X-large 3cm | 10/box |
| MAD062 | Adaptic Digit Non-Adherent Dressing Toe 2.8cm | 10/box |
| Packing- Use for packing wounds | | |
| Derma Science | | |
| DUP59145 | Iodoform Packing 0.25" x 5 yards | 12/box |
| DUP59245 | Iodoform Packing 0.5" x 5 yards | 12/box |
| DUP59345 | Iodoform Packing 1" x 5 yards | 12/box |
| Moisture Retention Hydrocolloid Use: To maintain wound moisture, these dsg do not absorb. Caution on elderly as it can cause skin stripping upon removal. | | |
| Convatec | | |
| 187955 | Duodern Extra Thin CGF 9.5cm x 9.7cm | 10/box |
| Wound Hydration Hydrogel Use: To add moisture to dry wounds, stop use when the wound begins to produce exudate Note: do not use in moist wounds | | |
| Smith & Nephew | | |
| 7308 | Intrasite Gel (8g Applicator) | 10/box |
| Antimicrobial Methylene Blue/Gentian Violet Use: To assist with removal of slough and eschar Note: Moisten dsg prior to application | | |
| Stevens | | |
| HB4414 | Hydrofera Blue Classic 10.2cm x 10.2cm | 10/box |
| HB6614 | Hydrofera Blue Classic 15.2cm x 15.2cm | 10/box |
| HBT0906 | Hydrofera Blue Classic 9mm Tunneling Dressing (1.2g), 2/package | 10/box |
| Polyester Mesh containing Silver Sulphate Use: Prevent wound trauma, low-mod exudating wounds | | |

| | | |
|--|--|--------|
| Urgo Medical | | |
| 509341 | UrgoTul Ag/Silver Contact Layer Dressing 10cm x 12.5cm | 10/box |
| Nanocrystalline silver mesh Use: Can be used for all type of wounds, burns and under NPWT | | |
| Smith & Nephew | | |
| 20101 | Acticoat 10cm x 10cm | 12/box |
| 66800545 | Acticoat Flex 7 Rope 2.5cm x 60cm | 10/box |
| 66800399 | Acticoat Flex 3 10cm x 10cm (3 day wear) | 10/box |
| 66800397 | Acticoat Flex 7 1cm x 12.5cm (7 day wear) | 10/box |
| Biochemical Modifier with Silver Use: to correct cellular imbalance in stalled wounds Note: do not apply on wound with eschar/slough | | |
| 3M | | |
| PS2028 | Promogran Prisma (Ag) 28cm ² | 10/box |
| PS2123 | Promogran Prisma (Ag) 123cm ² | 10/box |
| Silver Alginate: Use: Can help disrupt biofilm | | |
| Coloplast | | |
| COL3755 | Biatain Alginate Ag 5cm x 5cm | 30/box |
| COL3760 | Biatain Alginate 10cm x 10cm | 10/box |
| COL3765 | Biatain Alginate Ag 15cm x 15cm | 10/box |
| COL3780 | Biatain Alginate Ag 3cm x 44cm RIBBON | 10/box |
| Compression Use: Management of venous and mixed leg ulcer - Suggest Drs order and APBI before application | | |
| 3M | | |
| 3M2094 | Coban 2 Layer Short Stretch Compression 10cm x 2.7 m | 1/box |
| 3M2794 | Coban 2 LITE Layer Short Stretch Compression 10cm wide | 1/box |
| Zinc impregnated gauze Use: for chronic eczema and dermatitis, can be used with compression but not a substitute for compression. | | |
| Smith & Nephew | | |
| SN4948 | Viscopaste 7.5cm x 6m | 1/box |
| Contact Layer Paraffin Gauze Dressing low-adhesive paraffin tulle gras gauze dressing Use: on donor/graft sites (used by Plastics) | | |
| Smith & Nephew | | |
| SN7404 | Jelonet 10cm x 10cm | 10/box |
| Wound Hydration Hydrogel Use: To add moisture to dry wounds, stop use when the wound begins to produce exudate. Note: Do not use in moist wounds. | | |
| Smith & Nephew | | |
| 7308 | Intrasite Gel 8G - single use only | 10/box |

| Wound Cleansers and Other General Supplies | | |
|---|--|--------|
| Trudell | | |
| T168000 | Normal Saline NACL 0.9 STERILE (100mls/bottle) | 25/box |
| Sterasure | | |
| 33055-250DC | BIHOCL Pure Cleanse Hypochlorous Acid Wound Cleanser (250mls/bottle) | 24/box |
| Puritan Medical | | |
| HAR2515061PFDM | Wound Measurement Device STICK | 50/box |
| Derma Science | | |
| DUP59120 | Gauze Packing Strip, Plain Sterile, 0.25 in x 5 y | 12/box |
| DUP59220 | Gauze Packing Strip, Plain Sterile, 0.5 in x 5 y | 12/box |
| DUP59320 | Gauze Packing Strip, Plain Sterile, 1 in x 5 y | 12/box |