

LTC Provincial Wound Formulary June 2023

Effective June 2023, all wound care products can no longer be ordered through a pharmacy and should be ordered directly through the igility group purchasing wound care order site www.igilitygpo.ca

EXUDATE MANAGEMENT

Note: Do not use foam dressing on feet. Do not use with hydrogels.

Foam	Lite Silicone Foam Use: Scant amount of exudate. Indicated for up to 7 days or when exudate reaches within 0.5 cm from the edge of dressing.		
	66800835	Allevyn Gentle Border Lite 10 cm x 10 cm	10/Box
	66800840	Allevyn Gentle Border Lite 15 cm x 15 cm	10/Box
	Border Silicone Foam Use: Moderate-large amount of exudate and skin tears. Indicated for up to 7 days or when exudate reaches within 0.5 cm from the edge of the dressing.		
	66800270	Allevyn Gentle Border 10 cm x 10 cm	10/Box
	66800272	Allevyn Gentle Border 12.5 cm x 12.5 cm	10/Box
	66800959	Allevyn Gentle Border 17.1 cm x 17.9 cm	10/Box
	Border Silicone Foam Use: Moderate to high amount of exudate and skin tears. Quadrilobe/4 lobes. Indicated for up to 7 days or when exudate is visible in 3 out of 4 lobes.		
	66801068	Allevyn Life 12.9 cm x 12.9 cm	10/Box
	66801069	Allevyn Life 15.4 cm x 15.4 cm	10/Box
	66801070	Allevyn Life 21 cm x 21 cm	10/Box
	Post-Op Silicone Foam Use: Post-op for Low to Moderate amount of exudate or frequent dressing changes. Indicated to stay in place for up to 7 days or when exudate reaches within 0.5 cm of dressing edge.		
	66800900	Allevyn Gentle Border Post-op 10 cm x 20 cm	10/Box
	66800264	Allevyn Gentle Border Post-op 10 cm x 25 cm	10/Box
	Non-Border Silicone Foam Use: Moderate amount of exudate, requires secondary fixation. Indicated for up to 7 days or when exudate reaches within 1.5 cm of dressing edge.		
	66802129	Allevyn Non Border 10 cm x 10 cm	10/Box
66802131	Allevyn Non Border 15 cm x 15 cm	10/Box	
Non-Silicone Polyurethane Foam Use: For moderate amount of exudate. Indicated up to 7 days or when exudate reaches 0.5 cm of dressing edge.			
66000044	Allevyn Adhesive Classic 12.5 cm x 12.5 cm	10/Box	
66000045	Allevyn Adhesive Classic 17.5 cm x 17.5 cm	10/Box	
Alginate	Calcium Alginate Use: Packing of wounds and additional absorption to manage moderate - heavy exudate under cover dressing. Indicated up to 7 days or if exudate begins to pool in the wound or if cover dressing becomes saturated.		
	66000520	Algisite M 10 cm x 10 cm	10/ Box
	66000522	Algisite M 2 cm x 30 cm	5/Box

LTC Provincial Wound Formulary
June 2023

ANTIMICROBIAL

Provide broad action and can be highly effective in killing microorganisms. For wounds that improve, antimicrobial dressings should be continued for 14–21 days.

Antimicrobial

10% Povidone Iodine Non-Adherent Sheet Use: For wounds with minimal exudate, or to assist in the drying of wounds. Indicated up to 7 days or when fading of color occurs. Requires a cover dressing.

P01481	Inadine Dressing 5 cm x 5 cm	25/Box
P01512	Inadine Dressing 9.5 cm x 9.5 cm	25/Box

Cadexomer Iodine Use: For sloughy moist wounds and management of biofilm. Effective up to 3 days or when Iodosorb turns yellow/grey. Requires a cover dressing.

66060630	Iodosorb Ointment 10 gm tube	4/Box
66001286	Iodosorb Powder 3 gm sachet	7/Box

Silver Alginate Use: For moderate - heavy amount of exudate. Caution: Do not use on dry wounds or wounds with minimal wound exudate. Indicated up to 7 days or when cover dressing is saturated.

3780	Biatain Alginate Ag 3 cm x 44 cm	10/Box
3760	Biatain Alginate Ag 10 cm x 10 cm	10/Box
3765	Biatain Alginate Ag 15 cm x 15 cm	10/Box

Nanocrystalline Silver Mesh Use: Can be used for all types of wounds, burns, and under NPWT. Indicated for 3-7 days.

66800545	Acticoat Flex 7 Rope 2.5 cm x 60 cm	5/Box
66800399	Acticoat Flex 3 10 cm x 10 cm (3-day wear)	12/Box
66800397	Acticoat Flex 7 12.5 cm x 12.5 cm (7-day wear)	5/Box

Polyester Mesh containing Silver Sulphate Use: Prevent wound trauma, low-moderate amount of exudate. Indicated up to 7 days.

509341	UrgoTul Ag/Silver Contact Layer Dressing 10 cm x 12.5 cm	10/Box
--------	----------------------------------------------------------	--------

0.2% Polyhexamethylene Biguanide (PHMB) Use: Gauze dressing for wound packing, tunnels, open surgical wounds. Indicated up to 7 days or when cover dressing is saturated.

Z7831AMD	AMD PHMB Packing Strips 0.63 cm x 0.91.4 m (1/4 inch)	10/Box
Z7832AMD	AMD PHMB Packing Strip 1.27 cm x 0.91.4 m (1/2 inch)	10/Box
Z7833AMD	AMD PHMB Packing Strip 2.51 cm x 0.91.4 m (1 inch)	10/Box
Z3332	AMD Kerlix Rolls 11.4 cm x 3.7 m	1 roll

Methylene Blue/Gentian Violet Use: To assist with removal of slough/eschar.
Note: Pre-moisten dressing 5-10 mins prior to application. Requires cover dressing to keep this dressing moist.

HB2214	Hydrofera Blue Classic 5 cm x 5 cm	10/Box
HB4414	Hydrofera Blue Classic 10.2 cm x 10.2 cm	10/Box
HBT0906	Hydrofera Blue Classic 9 mm tunneling dressing	10/Box

Charcoal and Silver Use: For malodorous wounds; contains silver.
NOTE: Do not cut, product can be folded to size.

MAS105	Actisorb Ag 220 Charcoal 10.5cm x 10.5cm	10/Box
--------	------------------------------------------	--------

LTC Provincial Wound Formulary June 2023

Protective	Non-Adherent Cover Use: Cover and protect areas. Indicated up to 7 days. Note: Does not promote moist wound healing.		
	7135	Primapore Adhesive Non-Woven Wound Dressing 8.3 x 6 cm	50/Box
	66000317	Primapore Adhesive Non-Woven Wound Dressing 10 x 8 cm	20/Box
	66000318	Primapore Adhesive Non-Woven Wound Dressing 15 x 8 cm	20/Box
Super-absorbent	Super-Absorbent Dressing Use: For absorbing high amounts of exudate. Indicated for up to 7 days. Can use either side of the dressing, the dressing pad needs to overlap the edge of the wound. Do not cut the dressing.		
	PRD500-025	Kerramax Care Super-Absorbent Dressing 5 cm x 5 cm	10/Box
	PRD500-050	Kerramax Care Super-Absorbent Dressing 10 cm x 10 cm	10/Box
	PRD500-100	Kerramax Care Super-Absorbent Dressing 13.5 cm x 15.5 cm	10/Box
	PRD500-120	Kerramax Care Super-Absorbent Dressing 10 cm x 22 cm	10/Box
	PRD500-380-B10	Kerramax Care Super-Absorbent Dressing 20 cm x 30 cm	10/Box
Skin Protectant	Liquid Protective Film Use: For protection of intact or damaged skin from wound exudate. Provides an effective film barrier for up to 4 days when used alone or under non-adhesive dressings. Allow liquid to be fully dried before application of cover dressing.		
	66800709	No Sting Barrier Film Spray 28 ml (single person use)	1 bottle
	59420600	No Sting Barrier Film Wipe 1 ml	50/Box
	Adhesive Remover Use: to assist with the removal of adhesive dressing and tapes. Wipe off after use.		
	37436	Allkare wipe (adhesive remover)	50/Box
	Fabric Sheet Impregnated with Silver Use: Reduces skin-to-skin friction, wicks moisture from skin surface, manages odor, & provides antimicrobial action. Do not use as a cover for creams, pastes, powders, or ointments. Must have tail exposed. Remove sheet daily for skin care. Single piece can be re-used up to 5 days unless soiled.		
	7910	InterDry Silver Roll 25 cm x 366 cm	1 roll
	Antiseptic	Chlorhexidine Gluconate (CHG) 2% Solution Compounded with 70% Alcohol Use: To achieve and maintain dry, stable eschar e.g. arterial wounds. Apply with friction for 30 seconds and reapply every 3 days. Do not apply to open wounds.	
101.07		Soluprep antiseptic wipe – For larger area	100/Box
102.03		Soluprep antiseptic swab – For smaller areas	50/Box
Betadine Povidone Iodine 10 % Use: On dry or open eschar to maintain an intact covering of a wound and eschar. apply daily. Discard open bottle after 30 days.			
1034105		Povidone Iodine (Betadine) 115 ml	1 bottle
Wound Cleanser	BIHOCL Pure Cleanse Antimicrobial Cleansing Solution Use: A safe and non-cytotoxic, antimicrobial cleansing solution for wounds, skin, and burns. Use in place of saline to cleanse wounds with each dressing change.		
	33055-118DC	BIHOCL Pure Cleanse Hypochlorous Acid Wound Cleanser 118ml	1 bottle
	Normal Saline Use: For wound cleansing. Hold tip of the container 10-15 cm (4-6 inches) from the wound and squeeze solution over wound bed in a sweeping motion using all the solution. Discard the bottle after use.		
	T168000	100 ml Normal Saline- single use	25/case

LTC Provincial Wound Formulary

June 2023

Compression Wraps	Edema Wear Tubular compression stocking provides mild longitudinal compression (~10--20mm/Hg). Use: Stocking should be worn all day/all night and the stocking should be removed/ repositioned at least once a day. The stocking can be left in place for up to one week and washed in cold water and hang to dry. Can be reused up to 6 months.		
	CA600001	Edema Wear Small up to 45cm/ Navy (foot to knee). Each package contains 2 stockings.	1 package
	CA960001	Edema Wear Medium 45- 75cm/Yellow (foot to knee/groin). Each package contains 1 stocking.	1 each
	Compression Therapy: Non-Elastic/Short-Stretch Wrap Use: Can apply for up to 7 days given slippage does not occur and/or, if a wound is present, the wound exudate is managed. PLEASE NOTE: Only health care professionals who have successfully completed additional education for compression therapy (as per organizational policy) may apply compression wrap.		
	3M2094	Coban 2 Layer Short Stretch Compression 10cm x 2.7 m	1/Box
	3M2794	Coban 2 Lite Layer Short Stretch Compression 10cm wide	1/Box
Wound Contact Layer			
Non -Adherent Contact Layer	Petrolatum Impregnated Dressing Use: For prevention of wound bed trauma by decreasing adherence of the secondary dressing. Can be left in place for up to 7 days. Requires a cover dressing that covers the wound, plus at least 2 cm of peri-wound skin.		
	2012	Adaptic Non-Adhering Dressing 7.6 cm x 7.6 cm	50/Box
	2015	Adaptic Non-Adhering Dressing 7.6 cm x 20.3 cm	24/ Box
Wound Hydration			
Hydrogel	Intrasite Gel -Water based Amorphous Gel Use: Add/maintain moisture in wounds with necrotic tissue to enhance autolytic debridement. Cover with moisture retention dressing. Do not leave dressings in place for longer than 2 days.		
	7308	Intrasite Gel 8G - single use only	10/ Box
Other			
	DUP59901	Cotton Tip Applicators	100 Bag
	195-7144301	Hypafix Conformable Adhesive Retention Tape 5 cm x 10 m	1 roll
	195-7144302	Hypafix Conformable Adhesive Retention Tape 10 cm x 10 m	1 roll
	7143616	Easifix Cohesive 8cm x 20 m	1 roll

LTC Provincial Wound Formulary

June 2023

Prevention			
Requires pressure injury risk assessment and documentation by Nursing, Physio, or Occupational Therapy.			
Preventative	Offloading Heel Boots		
	MDT823330P	Heel Medix Advanced Petite 12.7 - 27.9 cm	1 boot
	MDT823330	Heel Medix Advanced Standard 27.9 - 48.2 cm	1 boot
	Positioning Wedge		
	MSC50023	Comfort Glide Positioning Wedge - 30° Foam Wedge	6 per case
	MSC50024	Comfort Glide Positioning Wedge (XL) - 30° Foam Wedge	2 per case
	Transfer Sheets		
	A 2-piece system. Need to order both the draw sheet and fitted sheet.		
	PTD-56	Swift Glide Draw Sheet (for mattresses 36" to 39" wide x 76" to 84" long)	1 sheet
	PTD/LS/S	Swift Glide Fitted Sheet	1 sheet
	PTD-56/BA/Y	Swift Glide Bariatric Draw Sheet (for mattresses 42" to 54" wide x 76" to 80" long)	1 sheet
	PTD/LS/S/BA/Y1	Swift Glide Bariatric Fitted Sheet	1 sheet



For more information contact group.purchasing@igility.ca



For wound care exceptions contact Woundcare@healthassociation.ns.ca